



Insert LP-Set for Piezo Z-Stage (Ludl) + CO₂-Cover LP


The Insert LP-Set consists of the frame and 4 base plates. The components marked with * can be used together with the set but are not included.




Insert LP-Set frame and base plate




Base plate for Petri dish 35 + POCmini





Base plate for Petri dish 60 + POC-R2



Base plate for Ibidi chamber




Base plate for Lab-Tek




Round and rectangular cell culture vessels


- Petri dishes 35mm and 60mm
- Cell Cultivation System POC-R2 and POCmini for open, closed culture and perfusion experiments
- Lab-Tek™ (Nunc)
- Chamber Slides (Falcon)
- normal slides
- IBIDI-Chambers




Base plate with Petri dish 35



Base plate with Petri dish 60



Base plate with Ibidi chamber




Base plate with Lab-Tek



CO₂-Cover LP transparent*



CO₂-Cover LP black*



Insert LP-Set with CO₂-Cover LP transparent*



Insert LP-Set with CO₂-Cover LP black*