



Insert PP-Set for Piezo Z-Stage (Prior) + CO₂-Covers PP




Insert PP-Set frame and base plate for Petri dish



Base plate for Lab-Tek




Base plate for Ibidi chamber




Frame and base plate with Petri dish 35



Base plate with Lab-Tek



Base plate with Ibidi chamber






Insert PP-Set with CO₂-Cover PP transparent*



Insert PP-Set with CO₂-Cover PP black*

The Insert PP-Set consists of the frame and 3 base plates.
The two CO₂-Covers marked with * can be used together with the set but are not included.



Round and rectangular cell culture vessels

- Petri dishes 35mm
- Lab-Tek™ (Nunc)
- Chamber Slides (Falcon)
- normal slides
- IBIDI-Chambers